



## RESmart GII BPAP System (T Series)

Bi-level positive airway pressure with Ti Min/Ti Max

RESmart GII BPAP Series provide non-invasive ventilation for patients with Obstructive Sleep Apnea (OSA) and respiratory insufficiency in hospital or at home. Ti Min/Ti Max ensure inspiration time and improve the treatment effect significantly.



## 3.5-inch colorful LCD Display

Real time pressure & flow wave display.

## SpO<sub>2</sub> Kit

Blood oxygen saturation and pulse rate display.

## Eco Smart Technology

Eco Smart heating system with innovative dual water chambers design, accurate water quantity control, with real time compensation. Eco Smart warm up, high efficiency, lower water bills.

## iCode QR


Compliance data is encoded by QR code. Search “BMCares” on Google Play or App Store.

## Ti Min/Ti max

Ensure inspiration time, improve the treatment effect significantly.

## Multiple data retrieval

SD Card / GPRS / Wi-Fi / iCode / iCode QR / iCode QR+



	T-20S	T-20A	T-20T	T-25S	T-25A	T-25T	T-30T
Mode	CPAP, S	CPAP, S, Auto S	CPAP, S, T, S/T	CPAP, S	CPAP, S, Auto S	CPAP, S, T, S/T	CPAP, S, T, S/T
Pressure		4 - 20 cm H <sub>2</sub> O			4 - 25 cm H <sub>2</sub> O		4 - 30 cm H <sub>2</sub> O
LCD Screen		2.4 inches				3.5 inches	
Weight		5.5 lb (2.5 kg) with humidifier			5.5 lb (2.5 kg) with humidifier		
Dimensions (L × W × H)		11.4" × 7.7" × 5.3" (290 × 196 × 134 mm) with humidifier			11.4" × 7.7" × 5.3" (290 × 196 × 134 mm) with humidifier		
SpO <sub>2</sub> Kit		Optional				Optional	
Data Storage		SD Card				SD Card	
Reslex		Off / 1 / 2 / 3				Off / 1 / 2 / 3	
Ramp P		4 cm H <sub>2</sub> O				4 cm H <sub>2</sub> O	
Max Ramp		0 - 60 minutes				0 - 60 minutes	
Ramp Time		0 - Max Ramp				0 - Max Ramp	
Humidity Level		off / 1 - 5				off / 1 - 5	
Ti Min		0.1 s - Ti Max				0.1 s - Ti Max	
Ti Max		0.3 s - 3.0 s				0.3 s - 3.0 s	

## BMC Medical Co., Ltd.

Add: Room 110 Tower A Fengyu Building, No. 115 Fucheng Road, Haidian, 100036 Beijing, PEOPLE'S REPUBLIC OF CHINA