



# Luna **Auto CPAP System**

E-20A-H-0/E-20AJ-H-0

With Luna, you will be fully supported on the path of sleep apnea treatment. E-20A-H-0 and E-20AJ-H-0 assembled with LCD screen and its 3-button control panel simplifies your life.



# Luna Auto CPAP System

# BMC

E-20A-H-O/E-20AJ-H-O

## Easy data access

iCode and SD card are two easy ways to transfer how you went with therapy. Those records can be viewed and evaluated through BMC+ iCode App or PC software.

# 01

# 02

## Smart pressure match(APAP)

APAP prevents upper airway collapsing with varied therapy pressure. Reslex respiratory pressure relief technology reduces the work of exhale so users remain comfortable and compliance.



## Other key features

Ramp | Mask test | Delay off | Multi-language

# 03

### BMC Medical Co., Ltd.

Add: 5/F Main Building, No.19 Gucheng Street West, Shijingshan, Beijing 100043, China

Tel: +86-10-51663880

Fax: +86-10-51663880 ext.810