



S01 5B Oxygen Concentrator



SO1 5B Oxygen Concentrator

Medical-standard oxygen supply

Provide oxygen with a concentration of 93%±3% for 24h.

Smart remote control

It allows you to set the device remotely (including backlight, operation time and oxygen flow rate), making oxygen inhalation easier.

Powerful kernel

Thanks to the stable and high-efficiency lithium molecular sieve, high adsorption efficiency, stable performance and sustainable high-concentration oxygen output can be achieved.



Multiple filtrations for clean breathing

Efficiently filter out impurities in the air to make the generated oxygen cleaner.

Convenient timed oxygen

Timed oxygen can be set up to 4 hours, and it will automatically shut down after the time is up.

More features are waiting for you to find out



Oxygen concentration displayed in real time



One-button start



High-definition large display



Visible humidifier



Universal caster with lock



Powerful cooling system



Replaceable filter cartridge

Specifications

SO1 5B Oxygen Concentrator			
Dimensions (L x W x H)	388 mm x 279 mm x 607 mm	Sound level	55 dB(A)
Weight	approx. 18 kg	Rating voltage	~ 230 V, 50 Hz
Rated flow	0.5 – 5.0 L/min	Rating power	360 VA
Output oxygen concentration	93% ± 3% (v/v)	Work mode	continuous flow, timed flow
Output pressure	50 ± 10 kPa		

BMC Medical Co., Ltd.

Room 110, Tower A Fengyu Building, No.115 Fucheng Road, Haidian, 100036 Beijing, PEOPLE'S REPUBLIC OF CHINA

Tel: +86-10-5166 3880 / Fax: +86-10-5166 3880 ext. 810 / en.bmc-medical.com