



Smart Connectivity • Towards the Future
EN-S9 Series Syringe Pump

Product Features:




Auto Loading

Fully automatic slider box, fast and precisely loading the syringe to ensure smoothly start of infusion.



Drug Library

Drug management system, which can save up to 5000 drugs, sortable with color, support DERS.




Flexible Clamp

The pipe clamp can be folded and rotated, compatible with horizontal and vertical poles.



Large Screen

7-inch true color capacitive touch screen, bring extraordinary display and operating experience.



Smart Connectivity

Department data synchronization management. Interact with workstation and central management system.



Relay Infusion

One-step operation, easy stacking, support relay infusion by IrDA.

Other Parameters:

Rate Range	VTBI	KVO Rate
0.01~2300mL/h	0.1~9999.99mL	0~10mL/h
Syringe Types	Occlusion Pressure	Syringe Empty Alarm
Support 1ml~60ml syringes	15 Levels, 50~1125mmHg	Support
Battery	Power Source	Waterproof
Lithium battery, 11.1V, 2600mAh	AC:100~240V, 50/60Hz DC:12~16V	IP44



Shenzhen ENMIND Technology Co., Ltd.

Add: 5th Floor, Block A, Defengsheng Building, No. 41 Dabao Road, Bao' an District 23, Shenzhen 518101, P.R. China

Web: www.enmind-tech.com
Email: sales@enmind-tech.com
Tel: 0086 755 2331 6007



Facebook

