



Smart Connectivity • Towards the Future
EN-V9 Series Infusion Pump

Product Features:




Electric control

Electric door lock and electric flow clamp. Comfortable and fast tube loading and unloading experience.



Drug Library

Drug management system, which can save up to 5000 drugs, sortable with color, supports DERS.




Flexible Clamp

The pipe clamp can be folded and rotated, compatible with horizontal and vertical poles.



Large Screen

7-inch true color capacitive touch screen, bring extraordinary display and operating experience.



Smart Connectivity

Department data synchronization management. Interact with workstation and central management system.



Relay Infusion

One-step operation, easy stacking, support relay infusion by IrDA.

Other Parameters:

Rate Range	VTBI	KVO Rate
0.1~2300mL/h	0.1~9999.99mL	0~10mL/h
Bubble Size	Cumulated Bubble	Occlusion Pressure
9 Levels, 10~800μL	50~1000μL/15min	15 Levels, 50~1125mmHg
Battery	Power Source	Waterproof
Lithium battery, 11.1V, 2600mAh	AC:100~240V, 50/60Hz DC:12~16V	IP44



Shenzhen ENMIND Technology Co., Ltd.

Add: 5th Floor, Block A, Defengsheng Building, No. 41 Dabao Road, Bao' an District 23, Shenzhen 518101, P.R. China

Web: www.enmind-tech.com
Email: sales@enmind-tech.com
Tel: 0086 755 2331 6007



Facebook

