



Butter fly 700/500 Led Shadowless Operation lamp



Product Characteristics



Brand Name	Topcare
Model Number	Butterfly 700/500
Bulb quantity	112/81
Bulb type	LED
Illuminance (Lux)	80,000-180,000/60,00-160,000
Average bulb life (hours)	= 80,000
Color temperature (K)	3700—5000
Diameter of light spot (mm)	160-350
Color rendering index	85-98
Lighting depth (mm)	=1200
Temperature rises at operator's head	< 1?
Mains voltage	AC100-240V 50/60HZ
Input power (W)	400