

PT LED 700/500

Led Shadowless Operation lamp



Product Characteristics

Technical parameters	LED700/500
Illumination (Lux)	40,000-180,000/30,000-160,000
Lamp bulb quantity	80/48
Bulb type	German OSRAM LED Bulb
Bulb average life	$\geq 80,000$ h
Colour temperature (K)	3700K ~ 5000K
Colour rendering index	85~98
Illumination depth(mm)	50-180cm
Facular diameter(mm)	160-280mm
Dimming Control	1%-100%
Temp Rise (doctor head)	$< 1^{\circ}\text{C}$
Input Power	AC100-240V 50/60HZ