

PT LED 700 Led Shadowless Operation lamp



Product Characteristics

Technical parameters	LED700
Illumination (Lux)	40,000-180,000Lux
Lamp bulb quantity	80
Bulb type	German OSRAM LED bulb
Bulb average life	=80, 000h
Colour temperature (K)	3700K ~ 5000K
Colour rendering index	85~98
Illumination depth(mm)	50-180cm
Facular diameter(mm)	160-280mm
Dimming Control	1%-100%
Bulb Power	80
Temp Rise (doctor head)	< 1℃
Input Power	AC100-240V 50/60HZ
Best Install Height(m)	2.7-3.1m