

TOPCARE

Handheld Pulse Oximeter

Brand: Topcare

Model: PT-100

Origin: China



Major Features

- ☒ Light weight for carrying & easy to use.
- ☒ Silicon rubber shell protection & stable bracket for table usage
- ☒ Using DB9 type connector compatible with Nellcor Spo2 Sensor.
- ☒ Support the oximeter probe for adult, neonate & infant.
- ☒ Big size 2.8-inch color TFT LCD display for Spot/PR/pulse bar graph/plethysmogram.
- ☒ Visual & sound prompt function.
- ☒ Adjust the parameters in friendly menu.
- ☒ Low Battery voltage indicator.
- ☒ Standard 4X AAA 1.5V Alkaline battery or rechargeable Li battery is available for power supply.

Topcare Technology

6F, Building F. Huaxinruiming Industrial Area, Dalang, Longhua Shenzhen, Guangdong, China.