



YK-6000A/A-T/B/B-T

Ultraviolet Phototherapy

Yonker is one of the biggest medical and healthcare product manufacturers in China.

Range of product: Patient monitor, Pulse oximeter, Blood pressure monitor,

Thermometer, Fetal Doppler, UV radiation, etc.

The whole production line is medical CE certified and mostly FDA certified.



TECHNICAL SPECIFICATIONS

Application Range

Applicable to clinical units for UV irradiation and treatment of

leucoderma, psoriasis, Pityriasis Rosea, eczema and other skin diseases.

Main technical indicators

Radiation intensity	UVA: 10-16 mw/cm2
	UVB: 3-16 mw/cm2
Working voltage	110V or 220V±22V, 50Hz±1Hz
Input power	≤40VA
Radiation source	Philips special UV lamp, with average life expectancy of 1000-1200 hours

Environment Requirements

Ambient temperature range	10-40 C
Relative humidity range	30%-75%
Atmospheric pressure range	700hPa-1060hPa
Transportation and storage temperature	-40-55 C
Transportation and storage humidity	≤95%

FEATURES

- With small size, lightweight, hand-held design, easy to operate and patients can take it along with themselves
- High radiation intensity, improving treatment effect
- Many display modes
- Humanized built-in timer design
- Comb for scalp psoriasis
- Suitable for home use
- Work distance: 3±0.5 cm

Dimension :

276mm*54mm*35mm

Weight (without battery):

650 g

Package including

- 1set main machine
- 1pcs power cord
- 1pcs manual
- 1pcs certificate

NOTE: The specifications are subject to changes without prior notice.

AUTHORIZED AGENT

Xuzhou Yongkang Electronic Technology Co., Ltd.

Address: 4F 8C Building Science and technology Area B Economic development district

Xuzhou, Jiangsu, China Post code: 221001

Tel: +86 516 87892766

Fax: +86 516 87892766-606

Email: office@yonker.cn

Website: www.yonker.cn

Yonker

YK-6000A/A-T/B/B-T

投诉

该页面通过 草料二维码 制作